# GLUMA® 2 Bond

|      |                                 |     |
|------|---------------------------------|-----|
| (GB) | Instructions for use .....      | 3   |
| (ES) | Instrucciones de uso .....      | 10  |
| (BR) | Instruções de uso .....         | 17  |
| (GR) | Οδηγίες χρήσης .....            | 24  |
| (CZ) | Návod k použití .....           | 31  |
| (HU) | Használati utasítás .....       | 38  |
| (LV) | Lietošanas instrukcija .....    | 45  |
| (LT) | Naudojimo instrukcija .....     | 52  |
| (PL) | Instrukcja obsługi .....        | 59  |
| (HR) | Upute za uporabu .....          | 66  |
| (RO) | Instrucțiuni de utilizare ..... | 73  |
| (TR) | Kullanma talimatı .....         | 80  |
| (RU) | Инструкция по применению .....  | 87  |
| (AE) | تعليمات الاستخدام .....         | 94  |
| (CN) | 使用说明 .....                      | 101 |
| (KR) | 사용 설명서 .....                    | 107 |
| (BG) | Инструкции за употреба .....    | 114 |
| (ID) | Petunjuk penggunaan .....       | 121 |
| (TH) | คำแนะนำในการใช้ .....           | 128 |

# GLUMA® 2 Bond

Instructions for use 

## Product description

GLUMA® 2 Bond is a light-cured single-component adhesive for use in adhesive restorative dentistry. GLUMA® 2 Bond has been developed for the adhesive bonding of plastic filling materials (e.g. composites, compomer, Polyglas®) to dental hard substances and laboratory-fabricated restorations (e.g. ceramics).

With GLUMA® 2 Bond, priming, bonding and desensitizing can be done in one step. GLUMA® 2 Bond is an ethanol based solution of light-activated, adhesive resins. Before using GLUMA® 2 Bond, the dental hard substance is conditioned with an etching gel. (Etch & Rinse-technique)

## Composition

GLUMA® 2 Bond contains:

- Methacrylate
- Ethanol
- Fillers
- Photoinitiators
- Glutaraldehyde

## Indications

- Adhesive fixing of direct composite, Polyglas® and compomer restorations
- Adhesive fixing of indirect laboratory-fabricated ceramic, Polyglas® and composite restorations (inlays, onlays, veneers and crowns)
- Treatment of hypersensitive tooth regions

Read the instructions for use carefully before use. Observe the safety instructions when using GLUMA® 2 Bond.

## Application

### 1. Adhesive fixing of direct composite, Polyglas® and compomer restorations

#### 1.1 Preparation

- The use of a rubber dam is recommended.
- Clean tooth with a toothpaste free from oil and fluoride.
- Prepare cavity as specified by the rules of adhesive filling treatment.
- Rinse prepared site with water and blow-dry.
- For deep cavities, we recommend protecting the dentine near the pulp with an appropriate base filling (e.g. calcium hydroxide preparation and glass ionomer cement base filling). Cover the dentine near the pulp only. Leave the remainder of the cavity free for adhesion of GLUMA® 2 Bond.

## 1.2 Conditioning

Prior to application of GLUMA® 2 Bond to the cavity preparation, pre-treatment of the tooth structure with an adequate etching gel is mandatory (e.g. Gluma Etch 35 Gel by Heraeus Kulzer).

## 1.3 Application of GLUMA® 2Bond

- Place GLUMA® 2 Bond in a Dappen glass, pick it up with an applicator tip or a soft disposable brush and apply it generously to the complete surface of the cavity.
- Allow GLUMA® 2 Bond to set for 15 seconds.
- GLUMA® 2 Bond contains ethanol as a solvent. Carefully evaporate solvent and residual moisture under a gentle air stream until no movement of liquid can be detected.

An excessively strong air flow at the start of blow-drying will dilute the bonding and may result in inadequate adhesion.

- The surface must be visibly shiny after application of GLUMA® 2 Bond and after evaporation of the solvent. Make certain that the complete surface of the cavity is fully covered. If the cavity does not appear to be completely shiny all over, apply GLUMA® 2 Bond a second time as described above. *Make sure that no excess remains in corners of the cavity.*

Polymerise GLUMA® 2Bond for 20 s with a halogen or LED polymerisation lamp. Use a Heraeus Kulzer Translux® polymerisation lamp or a polymerisation lamp with comparable intensity (min. 400–500 mW/cm<sup>2</sup>).

Insufficient lamp intensity will result in inadequate adhesion. Polymerisation lamps should be tested with reliable testing devices at regular intervals. The light emitting window should be placed as close as possible to the plastic surface during polymerisation.

- Then immediately apply and process the filling material in the cavity as directed by the manufacturer.

## **2. Adhesive fixing of indirect laboratory-fabricated ceramic, Polyglas® and composite restorations (inlays, onlays, veneers and crowns)**

### ***2.1 Preparation of the restoration***

Prepare restoration surface as directed by the manufacturer's instructions for use.

### ***2.2 Conditioning***

Proceed as described in 1.

### **2.3 Application of GLUMA® 2Bond**

Make sure that no excess material is left in the corners of the cavity and that the GLUMA® 2 Bond coating is carefully blow-dried before polymerisation.

Otherwise proceed as described in 1.

- Polymerise GLUMA® 2 Bond for 20 s with a halogen or LED polymerisation lamp (see instructions above in 1.3).
- Apply resin cement (self-curing or dual-curing resin cement) as directed by the manufacturer's instructions for use.

### **3. Treatment of hypersensitive tooth regions and cleaning the tooth**

- Carefully rinse the surface with water and blow-dry.

#### **3.1 Conditioning hypersensitive tooth regions**

Proceed as described in 1.

#### **3.2 Application of GLUMA® 2Bond**

Proceed as described in 1.

- Polymerise GLUMA® 2 Bond for 20 s with a halogen or LED polymerisation lamp (see instructions above in 1.3).
- Carefully remove the oxygen inhibition layer with an alcohol-soaked pellet.

If the GLUMA® 2Bond desensitisation effect is not sufficient, apply again as described above, cure and carefully remove the inhibition layer with an alcohol-soaked pellet.

### **Safety instructions**

- GLUMA® 2Bond contains ethanol. Ethanol is very flammable. Do not inhale vapour.
- GLUMA® 2Bond contains glutaraldehyde. If it is inhaled or swallowed, it may damage health. Seek medical attention immediately if swallowed.
- GLUMA® 2Bond contains methacrylate and glutaraldehyde. Contact with skin or mucous membranes may cause irritation or sensitisation. Prevent contact with skin, mucous membranes and eyes by taking precautions such as using a dental dam (patient), wearing rubber gloves (dentist) and wearing safety glasses (patient, dentist). Rinse eyes with water for several minutes after contact. Seek ophthalmological medical attention if pain persists. Wash skin with water and soap after contact.

### **Special instructions:**

- For use by dental personnel only.
- Keep out of reach of children
- Materials containing eugenol may affect the polymerisation of GLUMA® 2Bond
- Container must be at room temperature before use
- Do not use GLUMA® 2Bond after the date of expiry

**Side effects**

This product or one of its components may in particular cases cause hypersensitive reactions. If suspected, information on the ingredients can be obtained from the manufacturer. In the event of a skin reaction or an allergy, do not use this product.

**Storage**

Store at temperatures of 4–25 °C (40–77 °F). The product can be stored in refrigerator. Use the adhesive at room temperature only. Do not expose to direct sunlight. The shelf life will be reduced by unsuitable storage conditions and the function of the product may be affected. Close the container tightly immediately after use.

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Heraeus Kulzer GmbH
Heraeus Kulzer, LLC